



INNOVA
PRO



SCULPTOR pro M20

High density electromagnetic energy

SCULPTOR pro M20

Muscle stimulation by SBSW is used for body shaping in the following areas:

ABDOMEN, THIGHS, BUTTOCKS, LEGS and ARMS.

How does it work?

The applicators placed in the target area generate an electromagnetic field that induces muscle contractions. These contractions liberate free fatty acids that stimulate muscle tone to achieve muscle improvement.

A single treatment of 30 minutes is equal to 30.000 squat workout but without the effort of the workout.



INNOVA
PRO






INNOVA
PRO

Visible
results
without
effort

METHOD

The electromagnetic field focalizes through the skin to stimulate the muscles with continuous involuntary contractions adequate for body toning, also inducing apoptosis.



SCULPTOR ADVANTAGES

- ✓ 2 hands free applicators.
- ✓ Up to 30.000 contractions per treatment.
- ✓ No consumables.
- ✓ Results without sweating.





SBSW Technology, WITHOUT SPORT, WITHOUT EFFORT OR SURGERY

SBSW: Symmetrical Biphasic Sine Wave;

The electromagnetic field interacts with the motor neurons that trigger muscular contractions.

It's a non-invasive technology used to develop and strengthen muscles through electromagnetic stimulation.

Upgrade your muscles, tone your body!

It is the way to influence on those areas that are more resistant to achieve results with an adequate diet and minimum efforts.

You will obtain spectacular results.

The energy penetrates the deepest tissues, stimulating and training the muscles with high intensity with our specially designed preset programs in terms of different target areas and purposes.

A REAL REVOLUTION FOR BUTTOCK LIFTING WITHOUT SURGERY!

TECHNICAL SPECIFICATIONS

SCULPTOR pro M20

Product name:	SCULPTOR PRO M20
Field intensity:	3.0 Tesla \pm 20%
Indications:	Muscle development, body toning
Pulse width:	260 μ s \pm 20%
Source:	3500W
Field frequency	Symmetric Bifasic pulse, máx 50Hz
Mode:	Automatic and manual
Cooling system:	Water and forced air
Working temperature:	5°C to 40°C
Storage temperature:	-5°C to 50 °C
Relative humidity:	< 80%
Atmospheric pressure:	86-106 Kpa
Dimensions:	460x380x1130 mm



REF: PZ-100





INNOVA
PRO

TOGETHER
WE CREATE **BEAUTY**

CENTRAL HEADQUARTERS

Street 16
Warehouse 39-41-43
Industrial Park Base 2000
Lorquí - Murcia

968 603 256

 [innovapro_oficial](https://www.instagram.com/innovapro_oficial)  [InnovaPro](https://www.facebook.com/InnovaPro)

www.innovapro.es

info@innovapro.es

Technology
SBSW